



Prosthetic socket or calf
circumference [mm]

The height of the cover in [mm]





A [mm]

B [mm]

Foot manufacturer, type and size

For example : OB Talco 1C50 L29

Knee manufacturer and type

For example : OB 3R80



20 cm

35 cm

15 cm

30 cm

10 cm

25 cm

Adapter

- ☐ 30
☐ 34

Colour

- ☐ Skin Tone1 ☐ Grey
☐ Skin Tone3 ☐ White
☐ Skin Tone6 ☐ Black
☐ Graphite

Finish

- ☐ Matt varnish
(charged extra)

Tattoo (charged extra)

- ☐ Catalog number from
the our website.
e.g.: OB-1001

- ☐ Own

Check list :



Photo of the foot or shell.

Both the outer and inner edges of the shell
contour must be visible. The camera lens must be
perpendicular to the pyramid.



4 photos,

as shown in the images above. The camera must be placed perpendicular to the legs, in a low
position, and both legs must be centered on the photo. Please remember that all prosthetic components must
be completely visible, including the socket. Please remove shoes, socks and any object that covers the
prosthesis.






Photo of the prosthesis from the lateral.

Position the camera so the top of the foot shell looks like a straight line.



NOTES



- All our covers can be shortened from the bottom by several mm.
- In the factory state, the lower part of the cover is flat.
- In the case of the foot as shown in the picture, the lower part of the cover can be ground to the shape of the foot.
- If something is incomprehensible, look at our videos on the website in the menu - "VIDEOS"

