











About Us

We are a small business based in the South of England, with over 20 year's experience in the prosthetics and orthotics industry.

As a small business, we do not have the overheads of larger organisations, which enable us to be extremely competitive in the market and provide a product of exceptional quality.

We can't begin to understand the psychological effect of losing a limb, but what we do understand is why many affected individuals want to recreate their external appearance in as much detail as possible. We are totally passionate about every life like silicone cosmesis we create. Our high definition silicone Prostheses and Orthoses not only make a positive difference to everyday living but it can help to rebuild confidence.

The Silicone Centre Creating life like silicone

Why Are People Choosing The Silicone Centre?

THE BEST VALUE

We offer a choice of finishes to suit every budget. Our silicone experts will do everything they can to take care of all your individual requirements.

WIDE COLOUR RANGE

Our unique skin match system assures a close match to the skin tone using an infinite number of colours.

HIGHLY DETAILED

A life like finish can be achieved through the enhancement of the finer details such as freckles, veins, moles, skin pigmentation, surface hair. We can offer both silicone and acrylic nails, both compatible for nail polish, resulting in a natural appearance.

HANDMADE IN THE UK

All our silicone products are manufactured in-house at our clinic in the UK. In addition to being so realistic and life like, hand sculpted silicone products are more durable and stronger.



An Affordable Range of Finishes to Suit Every Budget. How Do They Compare?

UNIFORM

Unlike most off-the-shelf silicone products and orthotics, with their unnatural look and cheap plastic feel, our high-quality silicone prostheses and orthotics are individually created by hand, using different strengths of silicone to provide the perfect level of support and comfort.

- Affordable solution, great as a secondary prosthesis
- · Nails, split toes and skin textures
- Close colour match to skin tone, with over 100 plus skin colours to select from
- Range of silicone strengths and thickness to provide support and comfort
- Highly durable material
- · One colour

TONED

In addition to all the features of the Uniform, the Toned finish provides additional layers of colour to offer a more natural appearance.

- Perfect for those looking for a prosthesis that blends with their natural skin colouring
- Enhanced detailing, great solution for limbs that are exposed such as arms and legs and for those special occasions or warmer months
- Three colours



REALITY

The Reality offers the most realistic life like appearance to the most bespoke level of finish. Several days of extreme attention and passion goes into fabricating each silicone prosthesis to the highest detail.

- Perfect for those looking for the ultimate expression
- Extensive anatomical analysis for precise sculpting
- Natural creases, veins, freckles, moles and hair, are just some of the finer details available
- Most accurate colour matching system
- Translucent multi-layering for unparalleled detail
- One of a kind prosthesis fully customised to each individual
- · Multiple colours including acrylic nails





The Foundations of Creating a Lifelike Silicone Prosthesis

Our goal at The Silicone Centre is to work closely with the Prosthetist to fabricate custom made silicone prostheses to help individuals restore self-esteem, confidence and comfort.

There are four key elements which we believe are the fundamental steps in achieving a high-quality prosthesis. We provide our partners with a clear and simple guide to each stage of the process and our requirements, any deviations or omissions to the process could have can have an adverse effect to the time line and accuracy of your order.

1. An accurate impression is essential for the technician to replicate. The prosthesis will only be as accurate and as lifelike as the impression. We recommend that Alginate is used to take hand or foot impressions










2. Achieving an accurate colour match is just as important as the impression. The Silicone Centre will provide a set of over 100 colour swatches for remote colour matching, this must be undertaken in good natural lighting. The swatches are all individually numbered, this enables the technician to identify and recreate the exact colour back in the lab. Please make a note of the colour swatch number(s) and simply complete the Patient Information Order Form.



3. High Quality Photographs provide us with not only a colour reference for fabrication, but also serve as a visual for our technicians. Photographs of both the affected and sound side are required from a wide variety of angles and distances (please note fingernails and toenails should be unpainted when photographed). Please inspect the photographs with the patient, instruct the patient to select their favourite for colour and detail.













Contact Us

Distributed by OPC Health in Australia & New Zealand

T: +61 3 9681 9666 F: +61 3 9681 9366 E: sales@opchealth.com.au W: www.opchealth.com.au



Order forms available online at www.opchealth.com.au